

NEW TO THE ZOLL® TEMPERATURE MANAGEMENT PORTFOLIO

The IQool™ System

The **New Standard** in High-Quality Surface TTM

- IQool pads are skin-friendly, non-adhesive soft silicone, allowing for quick, easy skin checks
- Closed-loop system with internal temperature feedback and bacteriostatic coolant
- IQool pads are pre-filled and can be wiped clean to last the continuum of care
- Intuitive interface graphics
- Downloadable case data for insightful reviews
- Comprehensive staff training and 24/7 Clinical Help line



ADVANCED



NON-ADHESIVE



WIPE TO CLEAN



COST EFFECTIVE



**LASTS THE
CONTINUUM
OF CARE**



**THREE ZONES
OF TREATMENT**




**DISCONNECT
AND GO**

HEAD/NECK SIZES

	Cooling pads		Stabilizing insulation			
	Empty	Full	X-small	Small	Medium	Large
Weight (kg)	0.59 (x 2)	0.83 (x 2)	0.12	0.13	0.14	0.15
Length (cm)			28.5	30.0	31.0	32.0
Width (cm)			23.0	24.0	25.0	26.0

THIGH SIZES

	Cooling pads				Stabilizing insulation	
	Small/Medium		Large		Small/ Medium	Large
	Empty	Full	Empty	Full		
Weight (kg)	0.86 (x 2)	1.33 (x 2)	1.00 (x 2)	1.50 (x 2)	0.17 (x 2)	0.22 (x 2)
Length (cm)					80.0 (x 2)	94.0 (x 2)
Width (cm)					29.0 (x 2)	29.5 (x 2)



TORSO SIZES

	Cooling pads						Stabilizing insulation		
	Small		Medium		Large		Small	Medium	Large
	Empty	Full	Empty	Full	Empty	Full			
Weight (kg)	0.94 (x 2)	1.38 (x 2)	1.08 (x 2)	1.60 (x 2)	1.19 (x 2)	1.86 (x 2)	0.28	0.33	0.38
Length (cm)							110.0	117.0	141.0
Width (cm)							44.0	48.5	49.5

Technical support: (877) 225-7487

The IQool™ System should only be used as indicated and by trained healthcare professionals. Only IQool System Pads and accompanying neoprene insulation wraps should be used with the IQool System. Refer to IQool System Instructions For Use for warnings, cautions, contraindications, indications for use and additional information. Manufactured by BrainCool, distributed by ZOLL Medical. MCN TP 2402 0361

ZOLL®